KÜNG Biotech & Umwelt

Konsumstrasse 20, CH-3007 Bern
Tel. 031 357 53 73
valentin.kueng@kueng-biotech.ch
www.kueng-biotech.ch

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Anhang

Guidelines for the validation and application of alternative inactivation methods to heat inactivation using an autoclave

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01/07/2016/v21

1. Table of Contents

1.	Table of Contents	2
2.	Introduction.....	3
2.1	Context	3
2.2	Aim of this manual.....	3
2.3	Definition of important terms and fundamentals.....	3
3.	Steps during the validation process of alternative inactivation methods	7
3.1	Work flow for validation and routine application	7
3.2	Factors affecting the validation process.....	8
4.	Alternative inactivation methods.....	9
5.	Conclusion	12
5.1	Minimal requirements.....	12
5.2	Critical factors for process validation	12
6.	Examples for the validation of alternative inactivation methods	13
6.1	NaOH for the inactivation of bacteria	13
6.1.1	Inactivation resistance and killing kinetic.....	13
6.1.2	Interfering factors of sodium hydroxide inactivation.....	14
6.1.3	Testing a worst-case scenario using GS spores	15
6.2	Ultraviolet light (UVC_{254nm}) for the inactivation of GS spores	16
6.3	Combined approach: NaOH and UVC_{254nm}	17
6.4	Chlorine to inactivate viruses and phages	18
6.4.1	Inactivation resistance and killing kinetic.....	18
6.4.2	Interfering factors of chlorine inactivation	19
6.4.3	Testing of a worst case scenario using MS2 Bacteriophages.....	20
6.4.4	Testing of the worst-case scenario using Adenovirus in a cell culture system.....	21
6.5	Testing chemical indicators for UVC_{254nm} and chlorine inactivation	24
6.5.1	Potential for using erythrosine B as chemical indicator for UVC _{254nm} inactivation	24
6.5.2	Potential for using DPD as chemical indicator for chlorine inactivation.....	24
7.	Resources for process validation protocols.....	27
8.	FAQs.....	29
9.	References	30

2. Introduction

2.1 Context

When handling organisms in a laboratory, liquid and solid cultures containing living organisms are generated. According to the Containment Ordinance (CO¹), disposal of these cultures has to be performed in such a way that they do not represent a hazard to the population or the environment after being released from the 'contained' system into the sewage or the public waste disposal system. The preferred method of liquid and solid culture inactivation is heat inactivation using an autoclave. However, the CO allows the inactivation and disposal of solid cultures by external companies. Liquid cultures may be inactivated using alternative inactivation processes. Nevertheless, the commentaries on the CO state that the alternative process needs to be tested and validated by the user and must be approved by the federal authorities. Federal authorities are receiving increasing numbers of applications for approval, since autoclaves are expensive and alternative inactivation processes at a first glance can appear easier and more cost effective. However, there are many aspects and parameters to consider for a proper validation and the routine application of these methods.

2.2 Aim of this manual

In this manual, we aim to provide guidelines for the validation and (routine) application of alternative inactivation methods. This is essential, because various factors may determine the success of both the validation and the application of alternative inactivation methods. This manual will outline challenges during the validation process and highlight specific minimal requirements that should be fulfilled. The manual is intended to support laboratory users and the cantonal and federal authorities in deciding whether certain minimal requirements are met.

2.3 Definition of important terms and fundamentals

In order to understand the basics of inactivation, we will first provide an overview of important terms used in combination with the inactivation of organisms if not provided in the CO¹. The overall aim of inactivation is to obtain a sterile product, which can be disposed of safely.

Inactivation: Partial or full destruction of a given activity up to destruction of the microbiological system².

Sterile: State of being free from viable microorganisms². (See also section of sterility assurance level: SAL).

Sterilization: Process used to reach a sterile product².

Disinfection: Process of reducing the number of viable microorganisms by various physical and chemical methods. Disinfection usually does not lead to sterility³.

Decontamination: Removal of microbiological contamination or reduction to an acceptable level³.

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² Biotechnology – Laboratories for research, development and analysis. Guidance for handling, inactivation and testing of waste (EN12740:1999)

³ See footnote 2

Validation: Documented procedure for obtaining, recording and interpreting the results needed to show that a process will constantly yield a product complying with pre-determined specifications³.

Killing kinetic

A killing kinetic curve results from the connection of data points obtained by measuring the number of surviving organisms after a specific incubation time or radiation dose exposed to a specific inactivation method.

The example shows surviving *Geobacillus stearothermophilus* spores and vegetative *Staphylococcus aureus* after autoclaving for specific time intervals. In these cases, the killing kinetic curves are linear and therefore the D-value can be defined as the time to reduce the number of a specific microorganism by 1 log cycle i.e. $N_0/N_t = 10$. The D-value is specific to a temperature. When indicating D-values for an organism, the temperature at which the organism was tested must also be stated. This is usually done as a subscript of 'D' (Figure 1).


Figure 1: Killing kinetic of *G. stearothermophilus* spores and vegetative *S. aureus* inactivated by autoclaving at 121°C. Over time, samples are taken and cultured to monitor the log number of surviving organisms. The D-value is defined as the time it takes to reduce the number of organisms by 1 log cycle. The figure also shows that spores of *G. stearothermophilus* are more resistant to autoclaving than *S. aureus*, as indicated by a less steep killing kinetic curve.

Sterility assurance level (SAL)

Sterility means the absence of all viable microorganisms including viruses (Sharp 1995). At present, a sterility assurance level (SAL) of 10^{-6} is desired for sterilization procedures, i.e., a probability of not more than one viable microorganism in an amount of one million sterilized items of the final product. Since a SAL of 10^{-6} is difficult to measure experimentally, model situations need to be created, with the help of which conclusions can be drawn regarding the treatment conditions necessary to attain sterility meeting the SAL of 10^{-6} (see below). These are done using the reduction rate (D-value: Figure 1).

By extrapolating the reduction after sterilization, a theoretical overall performance of the procedure of 12 log increments (overkill conditions) is demanded to verify an SAL of 10^{-6} (see Figure 1). By comparison, other recommendations for thermal sterilization procedures demand only evidence that the difference between the initial contamination and the number of test organisms at the end of the process amount to more than six orders of magnitude. A practical proof of the required level of sterility assurance of 10^{-6} is not possible. For more details on sterility assurance levels, see von Woedtke and Kramer (2008).

Deviations from the sterility concept

Theoretically, when using non-thermic inactivation methods (anything other than the autoclave), killing kinetic curves are not linear and result in a so called 'tailing' (Lambert and Johnston 2000; Kramer and Assadian 2008). There are many hypotheses for this, but one

reasonable explanation for the tailing effect of killing kinetic curves with chemical disinfectants is 'quenching' i.e. the depletion of the disinfectant resulting in a reduced concentration (Johnston et al. 2000). In some cases, e.g. when excess chemical is used, nearly linear killing kinetic curves can be achieved. Nevertheless, it is important to note that a SAL of 10^{-6} cannot be guaranteed for the inactivation methods with non-linear killing kinetic curves.

Biological indicators

A biological indicator is a characterized preparation of a specific microorganism with a defined concentration that provides a defined and stable resistance to a specific inactivation process. If the inactivation resistance to a specific process is not known, it has to be demonstrated that it is more resistant than the organisms present in the unit of material to be sterilized. Biological indicators are usually viewed as a model for biological load to simulate the most resistant biological indicator organism⁴. Spore-forming microorganisms (e.g. *G. stearothermophilus*, *Bacillus atrophaeus* and *B. pumilus*) are widely used as biological indicators, because these microorganisms are highly resistant to inactivation compared to other microorganisms (Figure 2). The inactivation resistance hierarchy of microorganisms depicted below represents a rough guide and may also depend on the inactivation method.

There are three conditions for the use of an organism as biological indicator:

1. The inactivation resistance of the specific microorganism should be as high as possible towards the inactivation method. The inactivation of such highly resistant microorganisms encompasses all less resistant organisms, including most pathogens.
2. The specific microorganism ideally should not be pathogenic.
3. The specific microorganism needs to be culturable without great effort.


Figure 2: Order of inactivation resistance. This hierarchy considers broad classifications of microbial categories and it is considered a rough guide to general susceptibility patterns of microorganisms to sterilization. The hierarchy may depend on the inactivation method. Figure adapted from Favero and Bond (1991).

A biological indicator can be used to assist the process development and process validation steps of an alternative inactivation method (Figure 3) or to assess the performance of the inactivation equipment during routine monitoring.

⁴ For more information on requirements for biological indicators, see EN 11138: Sterilization of health care products - Biological indicators - Part 1: General requirements (ISO 11138-1:2006)

Types of biological indicators

There are three types of biological indicators. Each type incorporates a known species of a microorganism of a known inactivation resistance.

Type 1: Spores that are already added to a carrier (disk, strip, filter paper, glass, plastic or other material), packaged to maintain integrity and viability of the inoculated carrier. Carriers should not be degraded by the inactivation process. This type of biological indicator mimics a surface and can be used to validate surface sterilization processes.

Type 2: Spore suspensions consist of a known concentration and are directly inoculated on or into a unit of material to be sterilized. The stock spore suspension contains dormant (non-germinating) spores and is held in a non-nutritive liquid, e.g. 40% ethanol. Spore suspensions are probably the best option for validating the efficacy of chemical disinfectants.

Type 3: Self-contained indicators consist of a spore suspension in its own medium in ampules. They often contain a dye, which indicates positive or negative growth following incubation. This type is not suitable for testing chemical inactivation methods, because the chemical disinfectant cannot penetrate the ampules (no contact with the biological indicator).

G. stearothermophilus, *B. atrophaeus* and *B. pumilus* spores can be purchased from different companies, e.g. BAG Healthcare⁵.

Chemical indicators and physical sensors^{6,7}

Chemical indicators and physical sensors are used to monitor the inactivation process. Chemical indicators are designed to respond with a characteristic chemical change to one or more of the physical conditions during the inactivation process. In contrast, physical sensors are usually a piece of equipment (e.g. pH meter) indicating that a certain condition has been fulfilled. Consequently, the chemical indicator/physical sensor is process-specific and each inactivation method requires the use of different chemical indicators or physical sensors. For example, when using an acid for inactivation, it is not possible to use the same chemical indicator for acids as for bases (Table 1). Chemical indicators are designed to detect potential inactivation failures, but they do not test for sterility.

Chemical indicators have been classified by the Food and Drug Administration (FDA) into six classes depending on their action.

Class I: process indicators. They are intended to be used with individual items to demonstrate that the items have been exposed to the inactivation process and to distinguish between processed and unprocessed items. An example is autoclave tape.

Class II: dynamic air removal test (formerly called the Bowie-Dick test)

Class III: single parameter indicators (measure only one of the parameters of the inactivation process)

Class IV: multi-parameter indicators. These indicators reveal a change in one or more predefined process parameters based on a chemical or physical change resulting from exposure to a process. They provide much more information about the inactivation cycle than a Class III indicator.

Class V: chemical indicators. These are designed to react to all of the critical parameters over a specified range of inactivation cycles. Their performance has been correlated to the performance of the relevant biological indicator under the labeled conditions of use. This means that chemical integrators closely resemble biological indicators, but cannot be used to replace them.

⁵ <http://www.bag-healthcare.com/hygiene-monitoring/sterilisationsindikatoren/biologische-indikatoren>
http://www.spdceus.com/monitoring_sterilization_process.htm
http://www.pharmacopeia.cn/v29240/usp29nf24s0_c1035.html

⁶ For a good overview of chemical indicators, see: <http://www.sterislifesciences.com/Products/Process-Indicators/Chemical-Indicators.aspx>

⁷ For more information on requirements for chemical indicators, see EN 11140: Sterilization of health care products - Chemical indicators - Part 1: General requirements (ISO 11140-1:2005)

Class VI: Emulating Indicators. These are designed to confirm the presence or absence of specific time and temperature parameters during a cycle, and integrate all the critical parameters of the inactivation cycle (temperature, saturated steam and exposure time).

Interfering factors

Several physical and chemical factors influence the efficacy of an inactivation method. These include, for example external factors such as temperature, humidity or light exposure, and internal factors such as pH, water hardness or organic load (Cremieux 1986). Depending on the inactivation process, there is a specific set of interfering factors which applies. In general, the activity of inactivation methods increases as the temperature increases, but there is a point at which e.g. chemical inactivation agents can degrade, if heated too much. The temperature at which chemical inactivation agents degrade can be defined as an exclusion factor. Concerning pH, there are different mechanisms. For example, an increased pH improves the antimicrobial activity of agents containing aldehydes (Table 1), but decreases the activity of others such as hypochlorites. Relative humidity influences the activity of gaseous agents such as ethylene oxide. Water hardness reduces the reduction rate of some chemical inactivation processes, because divalent cations (e.g. Mg^{2+} and Ca^{2+}) interact with soap to form insoluble precipitates.

Organic load (proteins) such as serum, blood or fecal material may interfere with the activity of the inactivation process in at least two ways: First, a chemical reaction between the disinfectant and the organic load may result in a complex that is less germicidal or non-germicidal leaving less of the active agent to attack the microorganisms. Secondly, organic material may protect microorganisms from attack by acting as a physical barrier. For the same reason, another interfering factor can be biofilm formation. Biofilms are formed by bacteria and fungi and their presence make it more difficult disinfectants to penetrate the microbial matrix.

3. Steps during the validation process of alternative inactivation methods

3.1 Work flow for validation and routine application


Figure 3: Workflow and involved factors during process development, process validation and routine application of inactivation methods.

1st step: Process development

During the process development of a new process phase, it is important to assess the killing kinetic curve (Figure 1) of the inactivation method using the maximal biological load. Suitable biological indicators (see 1.2) should be used to demonstrate that the process is leading to the desired reduction rate (e.g. SAL of 10^{-6}). Suitable chemical indicators and/or physical sensors should be used to monitor the chemical/physical inactivation conditions. Different interfering factors (see 1.2, 5.1.2 and 5.4.2) need to be assessed. Process development only needs to be done once for an entirely new inactivation method. The result of the process development phase is a definition of the worst case scenario which is based on the worst possible combination of interfering factors. Therefore, all interfering factors first need to be identified during this step, although discovering all possible interfering factors is not always feasible. Using this worst case scenario, the validation process can be performed for new processes.

2nd step: Process validation

During the process validation, a reference load (=worst case scenario) is tested with biological indicators. The physical and chemical conditions are also monitored. Process validation needs to be done for new processes or first-time applications of a previously described (published) process. It is especially done when a new setup or a new instrument has been purchased.

3rd step: Routine monitoring

During routine monitoring of an inactivation method, a routine load is inactivated and biological and chemical indicators are used to assure functionality of the equipment and verify the physical and chemical conditions. Especially chemical indicators or physical sensors are used at every inactivation cycle. Inactivation of biological indicators needs to be verified at regular intervals.

3.2 Factors affecting the validation process

During all steps of process development and validation, the success of the inactivation method depends on whether **liquids or solids** are being inactivated. Therefore, different systems and different requirements are valid for solid and liquid wastes and these parameters need to be verified. The biological indicators and chemical indicators

Since microorganisms also differ in their inactivation resistance (Figure 2), it is necessary to keep in mind **what types of organisms** are being inactivated: **bacteria, spores, viruses or fungi**. The parameters used during process validation will differ depending on the organism. It is also important to keep in mind that published D-values are always organism and inactivation method specific and results are not transferable. Therefore, D-values must be determined for each species and at varying conditions considering interfering factors.

During all validation steps, the following treatment parameters may need to be monitored:

For most inactivation methods:

- Exposure time
- Volume of unit

For some inactivation methods:

- Concentration: important for chemical inactivation agents
- Temperature: important for chemical inactivation agents, but not so important for e.g. UV
- Pressure: important for gaseous agents
- Surface texture: important when sterilizing surfaces

4. Alternative inactivation methods

The selection of an appropriate alternative inactivation method may be determined by a number of criteria. These include:

1. Effectiveness of the alternative inactivation method
2. Applicability of the method to different organisms and media (solids or liquids)
3. Detoxification requirements and lack of toxic byproduct formation
4. Occupational and environmental hazards associated with the alternative inactivation method
5. Ease of handling and application
6. Operational costs

In this chapter, we provide a catalogue of different alternative inactivation methods considering the criteria mentioned above.

Table 1: Different inactivation methods and their mode of action, cost and advantage/disadvantage.

	Method	Active component	Spec-trum ⁸	Activity range	Mechanism	Known inter-fering factors	Combination options	costs	Comments Environmental/safety concerns
Physical approaches	Autoclave	Hot water vapor	Ba, Vi, Fu, Sp, *	100% rF	Denaturation	Not suitable for large volumes	not necessary	\$\$\$	High energy costs, not suitable for large volumes
	Ultraviolet Light (UV)	Energy-rich photons	Ba, Vi, Fu, Sp, *	250-300nm, 400J/m ²	Physical DNA damage	Turbidity (tailing effect)	all	\$\$	High energy costs, not effective against Listeria and Ameba
	Micro-/Ultra-Filtration	Sand filter, plastic filter membranes	Ba, Vi, Fu, Sp, *	NA	Mechanic separation	Turbidity, filters can clog	all	\$\$	Disposal of filter necessary, spectrum depends on pore size of filter. Spores are very small, only recommended when combined with other approaches, only small volumes
	Ultrasound	Low frequency Ultra-sound (20-35kHz)	Ba, Vi, Fu, Sp, *	500-5000 J/L	Cavitation, mechanic shearing forces on cells		all	\$\$	High energy costs, only for combined use

⁸ Abbreviations: Ba: bacteria, Vi: viruses, Sp: spores, Fu: fungi, *: reduced or no efficacy towards Prions

	Method	Active component	Spectrum ⁸	Activity range	Mechanism	Known interfering factors	Combination options	costs	Comments
									Environmental/safety concerns
	Ionized radiation	β^- and γ -radiation (^{60}Co , ^{137}Cs)	Ba, Vi, Fu, Sp, *	10^4 J/m^2	Protein and DNA damage		NA	\$\$\$	Radiation emitters are very expensive and process generates radioactive waste
Chemical approaches	Ozone	O_3 , H_2O_2	Ba, Vi, Fu, Sp, *	0.1mg/L	Oxidation	Organic compounds	UV	\$\$\$	Requires a lot of energy, generates toxic compounds, formation of hydroxyl-radicals in water, no storage possible, residual O_3 has to be removed
	Chlorine (halogens)	Cl_2 , NaClO , $\text{Ca}(\text{ClO})_2$, ClO_2	Ba, Vi, Fu, Sp, *	1-5%	Oxidation	Organic reactive compounds, light	filtration, Ultrasound, UV	\$	unstable, not biodegradable, corrosive
	Bases	NaOH	Ba, Vi, Fu, Sp, Prions	>1 Molar	Protein denaturation	Acids (neutralization)	heat, filtration, chem.	UV, \$	Not effective against Mycobacteria, low pollution potential, can be neutralized using acids, significant safety concern for user when handling product
	Acids	H_2SO_4 , $\text{C}_6\text{H}_8\text{O}_7$,	Ba, Vi, Fu, Sp, *	1 Molar	Protein denaturation	Bases (neutralization)	heat, filtration, chem.	UV, \$	Low pollution potential, can be neutralized using bases, significant safety concern for user when handling product
	Aldehyde	Formaldehyde, Glutaraldehyde, Glyoxal	Ba, Vi, Fu, Sp, *	0.5-5%	Protein denaturation	Development of resistance	heat, filtration, chem.	UV, \$	Stable, non-polluting, some compounds carcinogenic, irritating to mucosa
	Per-compounds	H_2O_2 , Peracetic acid (PAA), Potassium-peroxomonosulfate	Ba, Vi, Fu, Sp, *	0.02-35%	Oxidation	Unstable compounds, activity reduced by catalases and peroxidases	heat, filtration	UV, \$	Biodegradable, broad spectrum of activity
	Phenol-derivatives	P-Chlor-m-Kresol, p-Chlor-m-xlenol, o-Phenyl-phenol	Ba, Vi, Fu, *	0.1-5%	Protein denaturation	Organic compounds	Heat, filtration, chem	UV, \$	Safety hazard (neurotoxicity, carcinogenic), not biodegradable and toxic, not recommended.
	quaternary ammonium compounds	Benzalkonium chloride, Cetylpyridinium chloride, Didecyl-dimethyl-ammonium chloride	Ba, Vi, Fu, *	2-7%	Surface active, partially denaturing	Organic compounds, anionic soaps, CaCO_3 , iron	UV filtration, chem. Ultrasound	\$	Not biodegradable, "esterquats" are more easily biodegradable, especially the fatty acids

	Method	Active component	Spectrum ⁸	Activity range	Mechanism	Known interfering factors	Combination options	costs	Comments Environmental/safety concerns
Combined approaches	Advanced Oxidation	UV (182nm), Ozone, H ₂ O ₂ , catalysators	Ba, Vi, Fu, Sp, *		Cleavage of H ₂ O ₂ to hydroxyl-radicals	See individual approaches	Is combined		Superior to individual approaches, optimization of energy and chemical input, reduction of chlorine concentration possible
	Ultrasound chlorine	+ Ultrasound first or simultaneously	Ba, Vi, Sp		See individual approaches	See individual approaches	Is combined		Superior to individual approaches, optimization of energy and chemical input
	Ultrasound + UV	Ultrasound first for reduction of turbidity, UV penetration more effective	Ba, Vi, Fu, Sp		See individual approaches	See individual approaches	Is combined		Superior to individual approaches, optimization of energy and chemical input, tailing-effect reduced. UV dose or exposure time can be reduced
	Combination of different chemical approaches	Different combinations possible	Improved spectrum		See individual approaches	See individual approaches	Is combined		Additive, synergistic and antagonistic interactions possible, broader spectrum of activity towards different organisms
	UV + chemical approaches	UV first or simultaneously	Improved spectrum		See individual approaches	See individual approaches	Is combined		Superior to individual approaches, optimization of energy and chemical input
	Chemo-thermic disinfection	NaOH, heat and chlorine	Improved spectrum		See above (individual approaches)	See above (individual approaches)	Is combined		Chemical doses can be reduced by heat.

5. Conclusion

5.1 Minimal requirements

- Determination of a worst-case scenario for the unit to be inactivated (BI or suitable test organism with high inactivation resistance, interfering factors present in the unit).
- Process validation ([Figure 3](#), step 2) needs to be performed for all alternative inactivation methods and for new equipment using the worst-case scenario. Step 1 can be omitted for published killing kinetic curves for BI using known chemical indicators.
- If killing kinetic curves of an inactivation method are not linear, it has to be demonstrated that no organism of a 10^8 CFU/PFU suspension of a BI survives (this is just a suggestion...).
- If a waste contains several different microorganisms, inactivation data need to be obtained for each of the biological agents of concern using the identical conditions.

5.2 Critical factors for process validation

- Non-linear killing kinetic curves → SAL cannot be guaranteed
- Assessment of all possible interfering factors
- Availability of biological indicators for a given inactivation method
- Availability of chemical indicators for a given inactivation method

6. Examples for the validation of alternative inactivation methods

In order to outline in more detail the different steps of process validation and to highlight some challenges to be encountered when developing and validating a new alternative inactivation method, we provide experimental data for three alternative inactivation methods, two chemical methods:

- sodium hydroxide (NaOH, see Chapter 6.1),
- free chlorine as in sodium hypochlorite (Javel) or sodium dichloroisocyanurate (NaDCC, see Chapter 6.4),

and one physical method:

- Ultraviolet light (UVC_{254nm}, see Chapter 6.2),
-

and one combination of two methods (NaOH and UVC_{254nm}, see Chapter 6.3).

We outline some strategies for process development and validation of these three inactivation methods.

Nevertheless, if you plan to use one of the same inactivation methods, our data do not replace process development or validation in your laboratory. Interfering factors and test organisms may vary leading to a different worst-case scenario in your laboratory.

6.1 NaOH for the inactivation of bacteria

Sodium hydroxide (NaOH, base, see [Table 1](#)) was used to inactivate two different organisms: bacterial spores of *G. stearothermophilus* (GS) and vegetative cultures of *Staphylococcus aureus* (SA). *G. stearothermophilus* spores were chosen, because they are also used as biological indicators and are very difficult to inactivate by autoclaving.

6.1.1 Inactivation resistance and killing kinetic

In a first step of process development, it is essential to assess the killing kinetic and inactivation resistance of both the biological indicator and the test organism to demonstrate that the biological indicator is more resistant to the inactivation process than the test organism. In our example with NaOH as chemical disinfectant, we first assessed the inactivation resistance of GS spores and vegetative SA. Different concentrations of NaOH were incubated with GS spores or vegetative SA for different time intervals, i.e. 1, 10, 20, 30, 40, 50, 60 and 120 minutes and tested for the number of surviving organisms. By 120 minutes, the number of surviving organisms of GS and SA was below the detection limit for the chosen sodium hydroxide concentrations ([Figure 4](#)). After incubation, the chemical disinfectant has to be neutralized in order to eliminate residual antimicrobial activity⁹. Since sodium hydroxide is a base, it can be neutralized with the acid HCl-HEPES. Neutralization (pH=7) was verified using pH-indicator strips.

For the validation of chemical disinfectants, a catalogue of potentially suitable neutralization media can be found in the corresponding European Norms (EN, see Chapter 7).

⁹ **EN13727** – Chemical disinfectants and antiseptics – Quantitative suspension test for the evaluation of bactericidal activity in the medical area.

After neutralization, the inactivated bacterial suspension is serially diluted and plated on the appropriate medium and enumerated.

Using this approach, we found that GS spores and vegetative SA exhibited a different inactivation resistance to sodium hydroxide ([Figure 4](#)). GS spores were more resistant than vegetative SA. A 0.1M solution was sufficient to inactivate vegetative SA cultures, but the GS spores required at least a 4M solution and a longer incubation time to reach the desired SAL of 10^{-6} . Nevertheless, for the chosen organisms and within the tested range of sodium hydroxide concentrations, inactivation of GS and SA follows a linear killing kinetic. This means that an extrapolation is possible and the desired SAL of 10^{-6} is likely eventually reached. For GS this is after around 205 minutes and for SA around 130 minutes ([Figure 4](#)).


Figure 4: Linear killing kinetic.

Killing kinetic of the inactivation of *G. stearothermophilus* (GS, red) spores using 4M NaOH and vegetative *S. aureus* (SA, blue) using 0.1M NaOH. At the indicated time points a sample was removed and analysed for surviving bacterial entities. Plotting the number of surviving organisms in a logarithmic way, the killing kinetics are linear and can be extrapolated up to a SAL-value of 10^{-6} .

6.1.2 Interfering factors of sodium hydroxide inactivation

In a second step of process development, different interfering factors for the inactivation with NaOH were tested. The first one was inactivation **temperature** ([Figure 5](#)). Depending on where the inactivation is performed (at room temperature in a laboratory or in a colder room of 10°C or 15°C), other killing kinetic curves are monitored. In the case of 2M NaOH solution, we showed that decreasing temperatures reduce the efficacy of the disinfectant as chemical inactivation method against GS spores and much longer incubation times (>800 minutes compared to around 300 minutes) are required to achieve the desired SAL of 10^{-6} .

Organic load is another interfering factor. We simulated the presence of organic load in a test solution by the addition of bovine serum albumin (BSA) to a 2M NaOH solution ([Figure 6](#)). By the addition of 0.3% BSA, the incubation time is doubled (300 minutes compared to around 600 minutes), meaning that organic load also reduces the efficacy of sodium hydroxide inactivation.


Figure 5: Influence factor "temperature".

Inactivation of GS spores using 2M NaOH at different temperatures as interfering factors. At the indicated time points a sample was removed and analysed for surviving bacterial entities. At the lower temperature of 15°C (black), an increased incubation time of just under threefold was needed to reach the SAL-value of 10^{-6} as compared to the higher temperature of 20°C (red).


Figure 6: Influence factor "organic load".

Inactivation of GS spores using 2M NaOH (at 20°C) in the absence (red) or presence of organic load (3% BSA, blue). At the indicated time points a sample was removed and analysed for surviving bacterial entities. The reduction rate is reduced thus leading to a twofold increase in the required time to reach a SAL of 10^{-6} .

6.1.3 Testing a worst-case scenario using GS spores

A worst-case scenario most likely consists of a combination of interfering factors. Below, we provide an experimental example for the cumulative effect of interfering factors (Figure 7). We tested organic load combined with a reduced temperature (15°C). Combined interfering factors (0.3% BSA and reduced temperature) slow down the inactivation process even more than when only one interfering factor such as temperature is present.


Figure 7: Worst-case.

Effect of combined interfering factors (lower temperature, 15°C and organic load, 3% BSA, blue) on the inactivation of GS spores with 2M sodium hydroxide (NaOH). At the indicated time points a sample was removed and analysed for surviving bacterial entities. The reduction rate is reduced thus leading to a 3.5-fold increase in the required time to reach a SAL of 10^{-6} (as compared with the higher temperature of 20°C and without organic load, red). Interestingly, the effect of organic load was less pronounced on the inactivation if it was performed at lower temperature of 15°C (blue vs black; as compared with blue vs red in Figure 6).

It is important to note that the interfering factors presented here are only examples of the most obvious interfering factors. Not all interfering factors have been determined for the inactivation of GS spores by sodium hydroxide. Therefore, a final worst-case scenario applicable for all lab situations cannot be defined yet based on these data.

6.2 Ultraviolet light (UVC_{254nm}) for the inactivation of GS spores

UV light (UVC_{254nm}) kills by inducing DNA damage. This DNA damage is different in vegetative bacteria compared to bacterial spores, thus affecting inactivation resistance. When vegetative bacteria are exposed to UV light, thymine dimers are formed between adjacent thymine molecules in the same DNA strand ultimately leading to DNA damage. This DNA damage is very difficult to repair, but some bacteria are capable of it, e.g. *Deinococcus radiodurans*. In contrast, bacterial spores form a photoproduct called 5-thymine-5,6-dihydrothymine (TDHT) upon UV light exposure and resistance to UV light has been shown to be linked to TDHT removal (Munakata and Rupert 1972) or the presence of DNA coating small, acid-soluble spore proteins, for example present in *B. subtilis* (Setlow 1992). Therefore, bacterial spores are far more resistant to UV light compared to vegetative bacteria.

We exposed GS spores to different doses of UV light and the number of surviving GS spores was monitored. When combining the measurements of the number of surviving organisms and UV dose, there are two different extrapolation possibilities resulting in different reduction rates (Figure 8). Compared to sodium hydroxide inactivation, the UVC_{254nm} killing kinetic curve is not linear.


Figure 8: Ambiguous killing curve.

Inactivation of GS spores using different doses of UVC_{254nm}. At the indicated time points a sample was removed and analysed for surviving bacterial entities resulting in an ambiguous killing kinetic (a). There are two possibilities to extrapolate the killing kinetic curve (b; blue and red) resulting in different reduction rates.

6.3 Combined approach: NaOH and UVC_{254nm}

Another promising alternative inactivation approach is the combination of two different inactivation methods (Table 1). Especially radiation-based inactivation approaches are easily combined with e.g. chemical inactivation agents. Therefore, we also tested the two former inactivation methods, NaOH and UVC_{254nm}, together (Figure 9). The combination of a 0.1M NaOH solution with 0.2mW/cm² UVC_{254nm} resulted in a much steeper killing kinetic curve demonstrating that the combination of methods can be very effective.


Figure 9: Increase of efficacy by combined methods.

Killing kinetic curves of GS spores using UVC_{254nm} combined with NaOH (red) compared to both inactivation methods alone (black and blue). At the indicated time points a sample was removed and analysed for surviving bacterial entities. The resulting killing kinetic curve is much steeper for the combined method compared to UVC alone (blue), i.e. the combined treatment is more effective thus needing an incubation time 10-times shorter to reach the SAL-value of 10^{-6} .

6.4 Chlorine to inactivate viruses and bacteriophages

6.4.1 Inactivation resistance and killing kinetic

The third tested inactivation method was based on the action of free chlorine ions, as in sodium dichloroisocyanurate (NaDCC, Haztab, Guest Medical). Other widely used products that act via free chlorine are sodium hypochlorite (household bleach or Javel), chlorine dioxide and chloramines-T. Sodium hypochlorite was compared with NaDCC solutions and performed equally well for the inactivation of bacteriophages. Nevertheless, the advantage of NaDCC tablets over hypochlorite solutions is that they can be stored for a longer time without losing activity. However, once in solution, NaDCC loses free chlorine faster (Coates 1985). We confirmed this by measuring the amount of free available chlorine in NaDCC and sodium hypochlorite solutions. Chlorine products have a broad spectrum of antimicrobial activity, but the exact mechanism by which free chlorine destroys microorganisms has not been elucidated in detail. It may result from a number of factors including: oxidation of sulfhydryl enzymes and amino acids, ring chlorination of amino acids, loss of intracellular contents, inhibition of protein synthesis, oxidation of respiratory components, decreased adenosine triphosphate production, breaks in DNA and depressed DNA synthesis (Dychdala 2001).

The aim of these experiments was to test the efficacy of virus inactivation using a chlorine releasing product. Since experiments with tissue cultures and viruses (e.g. adenovirus) are very laborious and quantification relies on quantitative PCR, we chose the *E. coli* bacteriophage MS2 as a surrogate test organism. Bacteriophage MS2 was ideal, because it is easy to propagate, handle and count, and it exposes a high inactivation resistance to chemical disinfectants compared to other viruses and bacteriophages (Lehmann and Bansemir 1987). Nevertheless, bacterial spores were shown to be more resistant to chlorine inactivation than MS2 bacteriophage (Clevenger et al. 2007; Oie et al. 2011).

The efficacy of NaDCC was tested using a similar experimental setup compared to GS spore inactivation using sodium hydroxide. NaDCC solutions were incubated with MS2 bacteriophage suspensions for 5 to 120 minutes and neutralized using a sodium-thiosulphate based neutralization medium (Woolwine and Gerberding 1995). The number of plaque forming units (PFU) was enumerated on an *E. coli* soft agar overlay at different time intervals after addition of NaDCC. First, different concentrations of NaDCC were tested (Figure 10).


Figure 10: Non-linear or ambiguous killing kinetic.

Killing kinetic of MS2 bacteriophage in response to three different concentrations of sodium dichloroisocyanurate (NaDCC, 0.05% black, 0.01% red, 0.005% blue). Higher concentrations reduced the number of surviving organisms below the detection limit (2.5×10^5 PFU/ml). At the indicated time points a sample was removed and analysed for surviving phage entities resulting in a non-linear or ambiguous killing kinetic and no sufficient inactivation.

From these experiments, we conclude that chlorine inactivation of bacteriophage MS2 does not follow a linear killing kinetic, meaning a linear extrapolation is not possible. Therefore, reaching the desired SAL of 10^{-6} cannot be guaranteed for NaDCC, but to achieve satisfactory results, a much higher concentration of NaDCC needs to be used. This, however, can also be problematic especially considering the release of chlorinated products into the environment, which are highly corrosive and not biodegradable.

6.4.2 Interfering factors of chlorine inactivation

Four potential interfering factors of chlorine inactivation were tested. The first one was **organic load** – again simulated by the addition of BSA (Figure 11). Two different concentrations were tested (0.03% and 0.3%). We found that an addition of 0.3% BSA reduces the efficacy of chlorine for bacteriophage inactivation drastically resulting in nearly flat killing kinetic curves (=no killing). Therefore, the addition of 0.3% BSA represents a worst-case scenario in terms of organic load. As a reference, we determined the expected organic load in a typical laboratory waste originating from tissue culture using a Bradford Assay. We found that this typical laboratory waste contains an organic load of around 0.2-0.3% protein. Nevertheless, when validating chlorine releasing products, we recommend either to assume a worst-case scenario (=0.3% organic load) or to measure the organic load of a typical unit to be inactivated using a Bradford or Lowry Assay.


Figure 11: Influence factor "organic load".

Inactivation of MS2 bacteriophage with sodium dichloroisocyanurate (NaDCC 0.05%) in the absence (blue) or presence of organic load (0.3% red and 3% BSA black). At the indicated time points a sample was removed and analysed for surviving phage entities resulting in a non-linear or ambiguous killing kinetic and no sufficient inactivation. The presence of organic load reduced the efficacy of chlorine inactivation even more, especially at 3% BSA.

We then tested the inactivation **temperature** as second interfering factor. We defined that at our laboratory a room temperature of 18°C would represent lowest possible temperature. Comparing room temperatures of 18°C and 25°C, we did not observe any difference in inactivation efficacy of NaDCC and bacteriophage MS2 (data not shown). In a next step, we experimentally showed that a critical interfering factor for NaDCC is **age** (alias **shelf life**). We tested four differentially aged NaDCC solutions (fresh, 2, 7 and 11 week old solutions; Figure 12). While the 7-week old solution still slightly reduced the titer of an MS2 bacteriophage solution, an 11-week old solution did not show any inactivation. We therefore recommend using solutions which are not older than 3-4 weeks and to test for the amount of free available chlorine (FAC) in the solution before and after the experiment using e.g. the Dimethyl-4-phenylenediamine (DPD)-method (see 5.4.3).


Figure 12: Influence factor "age of solution/shelf life".

Inactivation of MS2 bacteriophage with 0.1% NaDCC of different age (fresh: black stippled; 2, 7 and 11 week old: black, red and blue solid, respectively). At the indicated time points a sample was removed and analysed for surviving phage entities. With the exception of the freshly prepared solution, all older NaDCC resulted in an insufficient inactivation exhibiting non-linear or ambiguous killing kinetics.

At the same time, we tested **light exposure** of these differentially aged solutions. We had stored an aliquot of each solution either exposed to light or in the dark (protected with foil). We found that in the containers used (50ml Falcon tubes), light exposure did not affect the disinfection efficacy of NaDCC (data not shown). Further, we tested whether light protected solutions of NaDCC still lost the same amount of free available chlorine over time and found that they did (data not shown). However, this might be different for other inactivation agents acting via free chlorine. For example, Rutala et al. (1998) found that storage of sodium hypochlorite in brown closed bottles significantly reduced free chlorine loss compared to solutions stored in clear containers. Another factor affecting free available chlorine loss over time is the concentration at which a solution is stored (Rutala et al. 1998). More concentrated solutions (1%) of sodium hypochlorite lose chlorine faster than more diluted ones (0.1%). Therefore, we recommend storing solutions at higher concentrations.

It is important to note that the interfering factors presented here are only examples of the most obvious interfering factors. Not all interfering factors have been determined for the inactivation of bacteriophages by free chlorine based products. Therefore, a final worst-case scenario applicable for all lab situations cannot be defined yet based on these data.

6.4.3 Testing of a worst case scenario using MS2 Bacteriophages

A worst-case scenario was defined as using a 6 week old NaDCC solution that contains organic load (0.3% BSA). We could clearly show that the efficiency of inactivation of MS2 bacteriophages for the old NaDCC solution were lower in the presence of organic load when compared to a freshly made NaDCC solution without organic load (Figure 13).


Figure 13: worst-case.

Inactivation of MS2 bacteriophage with 0.05% NaDCC freshly prepared without organic load ("best case") compared to 0.05% NaDCC 6 week old in the presence of 0.3% BSA. At the indicated time points a sample was removed and analysed for surviving phage entities. A more substantial loss in inactivation efficacy for the aged NaDCC solution with organic load (blue) compared to the "best case" (red) was observed as in the case of only organic loaded NaDCC (Figure 11, red curve). However, both cases NaDCC resulted in an insufficient inactivation exhibiting non-linear or ambiguous killing kinetics.

6.4.4 Testing of the worst-case scenario using Adenovirus in a cell culture system

The same worst-case scenario as described in 5.4.3 was used in an adenovirus (Ad5) infection experiment. HEK293 cells were infected with Ad5, which were differentially inactivated using fresh or NaDCC solutions aged in the presence of an organic load at different concentrations. The free chlorine activity was neutralized by cleaning the solution over a Sepharose column (Sephadex LH20; Figure 14). The increase in absolute numbers measured by qPCR of adenovirus genomes for each condition in the supernatant of the cell culture well indicates the presence of infectious adenovirus. For the experiment, we used different times of exposure [5 / 15 / 30 / 60 min], different NaDCC concentrations [1% / 0.1% / 0.01% / 0.001%], fresh vs. aged NaDCC solutions and presence or absence of organic load (0.3% BSA). In addition, several controls were added to test the functionality of the bioassay (Figure 15).

As shown in Figures 15, freshly made solutions are sufficient to completely inactivate the adenovirus at all tested concentrations, while NaDCC solutions aged in the presence of an organic load exhibited a significantly lower inactivation efficacy.


Figure 14: Neutralisation of the toxic effect of the inactivation solution using a Sepharose column.


Figure 15: Worst-case.

Inactivation of Ad5 with (a) freshly made NaDCC solutions of different concentrations (1%, 0.1%, 0.01%, 0.001%) without additional organic load in comparison with (b) NaDCC solutions of different concentrations (1%, 0.1%, 0.01%, 0.001%) aged 36 days in the presence of organic load (0.3% BSA). At the indicated time points a sample was removed and analysed for Ad5 infectivity. Green and violet bars indicate the total amount of Ad5 genomes 1 day (PI1) and 10 days after infection (PI10), respectively. In the worst-case scenario (b), except for the most concentrated solution (1% NaDCC), all concentrations failed to completely inactivate Ad5 exposed to the respective solution whereas in the "best case" (a) all NaDCC concentrations were sufficient to completely inactivate Ad5.


Figure 16: Controls.

(a) To test the linearity of the Ad5 reproduction in the infectivity assay (i.e. Ad5 genome copies in the supernatant) with or without prior Sepharose cleaning, different amounts of Ad5 (10^3 to 10^7) were added to infect HEK293 cells and the Ad5-genome copies in the cell culture supernatant were measured at days 1, 3 and 7 post infection.

(b) To test the functionality of the infectivity assay using infectious Ad5 in PBS (with or without additional organic load, 0.3% BSA) as well as heat-inactivated Ad5 samples or PBS itself were added to HEK293 and the cell culture supernatant was analysed for Ad5 genome copies.

6.5 Testing chemical indicators for UVC_{254nm} and chlorine inactivation

6.5.1 Potential for using erythrosine B as chemical indicator for UVC_{254nm} inactivation

UVC_{254nm} dosage is continuously measured with a radiometer. This radiometer represents an important physical sensor and should be used at every inactivation cycle using UVC_{254nm}. However, the measured dosage may not be the same in the entire load. Therefore, the use of additional indicators is recommended. Erythrosine B was tested as chemical indicator for UVC_{254nm} inactivation. Erythrosine B is traditionally used in the food industry as food coloring, and its absorbance peaks at 525nm (Figure 17 a). This absorbance decreases largely following an exponential function to around 700 mJ/cm² (Figure 17 b; Putt et al. 2012). With higher dosages than 700 mJ/cm², the decrease in percentage of remaining absorbance is not linear anymore. During inactivation, cuvettes containing erythrosine B are placed within the unit. After inactivation with a specific UV dosage, the absorbance at 525 nm is measured using absorptimetry. The color change from red to colorless can also be seen by the naked eye (Figure 17 c).


Figure 17: (a) Absorbance spectrum of erythrosine B reaching a maximum at 525nm. (b) Percent of remaining absorbance of erythrosine B exposed to increasing dosages of UVC_{254nm}. The absorbance drops exponentially until to approximately 700 mJ/cm². (c) Color change of erythrosine B upon exposure to UVC_{254nm}.

A major advantage of this system is that different cuvettes can be placed at different positions within a unit yielding values for the emitted UVC_{254nm} dosage at different locations in varying distances to the UVC_{254nm} source. In conclusion, we found that erythrosine B has potential to be used as chemical indicator for UVC_{254nm} inactivation, but still needs to be verified further on a larger scale more closely resembling the laboratory setting (e.g. testing with larger volumes or liquids and solids).

6.5.2 Potential for using DPD as chemical indicator for chlorine inactivation

We evaluated the potential of DPD (N, N-diethyl-p-phenylenediamine sulfate) as chemical indicator to measure the amount of free available chlorine (FAC) in NaDCC or sodium hypochlorite solutions. The amount of FAC is the amount of chlorine which is not bound to nitrogenous compounds (=organic load) and which is able to react with microorganisms. FAC is the determining factor for inactivation. Chlorine reacts with DPD (Oxycon DPD, Figure 18 a) in the presence of a suitable buffer (Oxycon Start, Figure 18 a) to produce a pink color (Figure 18 b, c, d). Depending on the concentration of FAC in a solution, the pink color is more intense, which can then be measured at 520 nm, using an LED-based spectrophotometer (Figure 18 e) yielding the amount of FAC in mg/L. Besides the amount of FAC, the DPD method can also be used to determine the amounts of combined chlorine (i.e. the amount of chlorine bound to a nitrogenous compound) and total chlorine using the Oxycon 2 buffer.


Figure 18: (a) Reagents used for chlorine measurements – Oxycon Start: buffer, Oxycon DPD: DPD for FAC measurements, Oxycon 2: for total and bound chlorine measurements (b) chlorine solution containing small amounts of FAC (<1000 mg/L in a 1% solution) (c) chlorine solution containing intermediate amounts of FAC (4000 mg/L) (d) chlorine solution containing high amounts of FAC (>6000 mg/L) (e) spectrophotometer for FAC measurements.

It has been shown previously that chlorine solutions cannot be stored indefinitely (Rutala et al. 1998), but some compounds lose FAC faster than others. In order to address this question in our NaDCC (Haztab) solutions, we monitored the concentration of FAC over time (for 36 days). We demonstrate that NaDCC solutions lose FAC over time and interestingly, the decline is nearly linear (Figure 19). We additionally monitored the concentration of FAC in NaDCC solutions in the presence of an organic load (0.3% BSA). This simulates the scenario where an inactivation solution is added to a container prior to the input of contaminated liquid waste. With the contaminated liquid waste, the organic load is introduced, and therefore the chlorine solution ages with the organic load. When aged together, a much faster (4x) decrease in FAC can be observed compared to just aging without organic load (Figure 19).


Figure 19: Process monitoring.

Free available chlorine in NaDCC (HazTab) solutions over time. The red line represents aging of a NaDCC solution without organic load added. The blue line indicates aging of a NaDCC solution in the presence of organic load (0.3% BSA).

NaDCC used here is only one example of chlorine solutions, normal household bleach (aka Javel, sodium hypochlorite) is frequently used in laboratories for the inactivation of contaminated liquid waste. To compare NaDCC with sodium hypochlorite solutions, we evaluated the amount of FAC and their inactivation efficacy. We measured equal starting concentrations of FAC in NaDCC and sodium hypochlorite solutions, but did not observe a decay of FAC (aged sodium hypochlorite solutions still had the same concentration of FAC after 36 days, data not shown). Nevertheless, they differed slightly in their inactivation efficacy (Figure 20).


Figure 20: Inactivation efficiency of a NaDCC solution and sodium hypochlorite solution (Axxel Javel: purchased from pharmacy) both diluted to a working concentration of 0.01%.

In conclusion, the DPD method can be used to quantify the amount of FAC in a solution. It is recommended to measure the amount of FAC in a solution prior to inactivation. When diluting chlorine solutions, the final concentration can be adjusted with the measurements obtained for FAC. Nevertheless, a disadvantage of the DPD method is that swimming pools contain far less FAC and thus, chlorine solutions need to be diluted. Depending on the amount of FAC, this can be up to 10^{-5} . Creating a dilution series can introduce errors in FAC measurements and subsequent calculation of FAC amounts.

7. Resources for process validation protocols

There are European Norms (EN), i.e. certified standard procedures which can assist during the validation process. They provide detailed protocols and represent helpful guidelines for the validation of different chemical disinfection methods and different organisms. They can be purchased for around CHF 150.00 from the publisher "Beuth" (www.beuth.de).

List of potentially useful European norms:

EN12740 – Biotechnology – Laboratories for research, development and analysis - Guidance for handling, inactivating and testing of waste (EN 12740:1999)

EN 11138 – Sterilization of health care products – Biological indicators – Part1: General requirements (ISO 11138-1:2006)

EN 11140 – Sterilization of health care products – Chemical indicators – Part 1: General requirements (ISO 11140-1:2005)

EN 1040 – Chemical disinfectants and antiseptics – Quantitative suspension test for the evaluation of basic bactericidal activity of chemical disinfectants and antiseptics-Test method and requirements (phase1)

EN 1275 – Chemical disinfectants and antiseptics – Quantitative suspension test for the evaluation of basic fungicidal or basic yeasticidal activity of chemical disinfectants and antiseptics – Test method and requirements (phase 1)

EN 1650 – Chemical disinfectants and antiseptics – Quantitative suspension test for the evaluation of fungicidal activity of chemical disinfectants and antiseptics used in food, industrial, domestic and institutional areas – Test method and requirements (phase 2, step 1)

EN 1656 – Chemical disinfectants and antiseptics – Quantitative suspension test for the evaluation of bactericidal activity of chemical disinfectants and antiseptics used in veterinary field-Test method and requirements (phase 2, step 1)

EN 1657 – Chemical disinfectants and antiseptics – Quantitative suspension test for the evaluation of fungicidal or yeasticidal activity of chemical disinfectants and antiseptics used in veterinary area – Test method and requirements (phase 2, step 1)

EN13610 – Chemical disinfectants – Quantitative suspension test for the evaluation of virucidal activity against bacteriophages of chemical disinfectants used in food and industrial areas (phase 2/step 1)¹⁰

EN 13624 – Chemical disinfectants and antiseptics – Quantitative suspension test for the evaluation of fungicidal activity of chemical disinfectants for instruments used in medical area – Test method requirements (phase 2, step 1)

EN 13697 – Chemical disinfectants and antiseptics – Quantitative non-porous surface test for the evaluation of bactericidal and/or fungicidal activity of chemical disinfectants used in food and industrial, domestic and institutional areas – Test method and requirements without mechanical action(phase 2,step 2)

EN 13704 – Chemical disinfectants – Quantitative suspension test for the evaluation of sporidicidal activity of chemical disinfectants used in food, industrial, domestic and institutional areas – Test method and requirements (phase 2, step 1)

¹⁰ This EN was used as a guideline for the validation of Javelle as chemical disinfectant.

EN13727 – Chemical disinfectants and antiseptics – Quantitative suspension test for the evaluation of bactericidal activity in the medical area (phase 2/step 1)¹¹

EN 14204 – Chemical disinfectants and antiseptics – Quantitative suspension test for the evaluation of mycobactericidal activity of chemical disinfectants and antiseptics used in veterinary area – Test method and requirements (phase 2, step 1)

EN 14347 – Chemical disinfectants and antiseptics – Basic sporicidal activity - Test method and requirements (phase 1)

EN 14348 – Chemical disinfectants and antiseptics – Quantitative suspension test for the evaluation of mycobactericidal activity of chemical disinfectants and antiseptics in the medical area including instrument disinfectants – Test method and requirements (phase 2, step 1)

EN14476 – Chemical disinfectants and antiseptics – Quantitative suspension test for the evaluation of virucidal activity in the medical area – Test methods and requirements (phase 2/step 1)

EN 14562 – Chemical disinfectants and antiseptics – Quantitative carrier test for the evaluation of fungicidal or yeasticidal activity for instruments used in medical area – Test method requirements (phase 2, step 2)

EN 14885 – Chemical disinfectants and antiseptics – Application of European Standards for chemical disinfectants and antiseptics

¹¹ This EN was used as a guideline for the validation of NaOH as chemical disinfectant.

8. FAQs

Q: If data concerning killing kinetic and reduction rate of my desired inactivation method with my desired organism is already published, do I have to redo these experiments?

A: As long as the same organisms, the same experimental conditions and the same interfering factors were present in the published dataset, the process development step can be skipped and directly proceeded to the process validation step to demonstrate validity at your laboratory.

Q: What do I have to do, if the killing kinetic of my desired inactivation method is not linear?

A: In this case, a much higher concentration of disinfectant and/or a much longer incubation time need to be followed in order to assure sterility. A SAL of 10^{-6} can never be guaranteed in this case, but absence of viable microorganisms needs to be demonstrated using biological indicators and validated assuming a worst-case scenario.

Q: How do I determine the organic load in my liquid waste, which I want to inactivate?

A: A Bradford Assay can be done to determine the protein quantity in your waste. A typical tissue culture waste contains around 0.2-0.3% protein.

Q: What is important when choosing biological indicators?

A: First you have to show that the biological indicator is more resistant than the organisms you want to inactivate. We would recommend using bacterial spores and comparing these results with the organisms most likely present in your inactivation material. This will show whether you can use commercially available biological indicators.

Q: What do I do, if there are no biological indicators for my desired inactivation method?

A: You cannot use this method.

Q: What do I do, if there are no chemical indicators for my desired inactivation method?

A: This is definitely not an ideal situation. Either you have to show that it is a validated product stored under the correct condition (always containing a concentration of XYZ). In addition, you have to use biological indicators at every inactivation cycle.

Q: How long do process development and process validation (as shown in Figure 2) take?

A: You should anticipate 6 months of experimental work to assess killing kinetic, test the interfering factors and test the chemical and biological indicators. Make sure to start with the latter two, so you do not waste time with the first two and then figure out that there are no chemical/biological indicators for your desired inactivation method.

9. References

- Clevenger, T., et al. (2007). "Comparison of the inactivation of *Bacillus subtilis* spores and MS2 bacteriophage by MIOX, ClorTec and hypochlorite." J Appl Microbiol **103**(6): 2285-2290.
- Coates, D. (1985). "A comparison of sodium hypochlorite and sodium dichloroisocyanurate products." J Hosp Infect **6**(1): 31-40.
- Cremieux, A. (1986). "Factors affecting the bactericidal action of disinfectants. Implications for selection of resistant strains." Drugs Exp Clin Res **12**(11): 899-903.
- Dychdala, G. R. (2001). Chlorine and chlorine compounds. Disinfection, sterilization and antisepsis in healthcare. S. S. Block. Philadelphia, Williams & Wilkins: 135-157.
- Favero, M. S. and W. W. Bond (1991). Chemical disinfection of medical surgical material. Philadelphia, PA, Lea & Febiger.
- Johnston, M. D., et al. (2000). "One explanation for the variability of the bacterial suspension test." J Appl Microbiol **88**(2): 237-242.
- Kramer, A. and O. Assadian (2008). Wallhäussers Praxis der Sterilisation, Desinfektion, Antiseptik und Konservierung. Stuttgart, Georg Thieme Verlag.
- Lambert, R. J. and M. D. Johnston (2000). "Disinfection kinetics: a new hypothesis and model for the tailing of log-survivor/time curves." J Appl Microbiol **88**(5): 907-913.
- Lehmann, L. H. and K. P. Bansemir (1987). "Bakteriophagen als Testviren für die Desinfektionsmittel-Prüfung - 1. Mitteilung: Quantitativer Suspensionstest mit Desinfektionswirkstoffen." Hygiene + Medizin **12**: 9-12.
- Munakata, N. and C. S. Rupert (1972). "Genetically controlled removal of "spore photoproduct" from deoxyribonucleic acid of ultraviolet-irradiated *Bacillus subtilis* spores." J Bacteriol **111**(1): 192-198.
- Oie, S., et al. (2011). "Disinfection methods for spores of *Bacillus atrophaeus*, *B. anthracis*, *Clostridium tetani*, *C. botulinum* and *C. difficile*." Biol Pharm Bull **34**(8): 1325-1329.
- Putt, K. S., et al. (2012). "The use of chromophore and fluorophore degradation to quantitate UV dose: FD&C dyes as chemical indicators for UV sterilization." J Microbiol Methods **91**(2): 215-221.
- Rutala, W. A., et al. (1998). "Stability and bactericidal activity of chlorine solutions." Infect Control Hosp Epidemiol **19**(5): 323-327.
- Setlow, P. (1992). "I will survive: protecting and repairing spore DNA." J Bacteriol **174**(9): 2737-2741.
- Sharp, J. (1995). "What do we mean by "sterility"?" PDA J Pharm Sci Technol **49**(2): 90-92.
- von Woedtke, T. and A. Kramer (2008). "The limits of sterility assurance." GMS Krankenhhyg Interdiszip **3**(3).
- Woolwine, J. D. and J. L. Gerberding (1995). "Effect of testing method on apparent activities of antiviral disinfectants and antiseptics." Antimicrob Agents Chemother **39**(4): 921-923.

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⁵⁰ Desinfektionsmittelliste des Verbund für Angewandte Hygiene (VAH), Stand: 1. April 2014; ©mhp-Verlag 2015; http://www.vah-online.de/uploads/PDF/vorwort_deutsch_mhp.pdf